

Protocol for Herring AOI images

August 2017

Camera: Leica DFC320 and Stereomicroscope: Leica MZ6

1. Settings for light box:

Images are taken under the circular reflective light.

Top knap 3.0 Bottom knap 3.0

2. Settings in LAS V4.5:

Setup tab

Microscope Configuration should NOT be changed!!!

Acquire tab:

MZ6 Be sure the microscope magnification is 2 and that the corresponding “ZoomDrive” is set at 2.

Acquire tab:

Camera

Input options

Current Configuration: “Sild_Circ_Aug 2017”

The saved settings are:

| Exposure Adjust | | Image Formats | |
|-----------------|-------|-----------------|--------------------------|
| Exposure | 414.5 | Captured format | 2088 X1550 Full Frame HQ |
| Gain | 1.0 | Live format | 2088 X 1550 Full Frame |
| Saturation | 0.00 | Image type | Colour |
| Gamma | 0.80 | | |

Acquire tab:

Camera

Calibration settings

Type: **Measured**

Configuration: **Sild_Aug_2017**

Actual length of line shown on image: 1

3. Set Capture location:

Images should be saved here: <O:\Sildeotolithbilleder til AOI>

Folders created by survey/collection type and then year

Use Navigator to find the folder where images will be saved

Press “Set Capture Location” icon


Taking images

1. Calibration Image


A calibration image needs to be taken each time a set of images are taken and saved in the same folder. If you change survey/collection type you have to take a new calibration image.

- Using the 5cm measuring stick focus the image until it is sharp and the marks appear clearly in the image viewer (see image below)

- The calibration image is named using “cal_0_5x1x2_0x1_date-month-year (using the format 12-07-17)”



- To take image: press “Acquire Image “



- In the “Basic Annotation” box place a check mark in “Show”
- Click “Recall” and then “Yes”
- Click “Merge All” and the “Replace”

2. Otolith images

- Images are taken of a pair of otoliths immersed in alcohol and on a black background. Be sure that the otoliths are not dirty (especially at the edges).
- It is VERY important that the following format is used for naming the images:
e.g. J0225_00001_0_5x1x2_0x1_ (following the format Jnr No_Fish No_Magnification)

You are then ready to take a set of images from a sample